SCIENTIFICALLY
PROVEN

pulse^{roll}

PERCUSSION MASSAGE GUN MANUAL

OFFICIAL SUPPLIERS

BRITISH
WEIGHT
LIFTING

DEAR CUSTOMER

I would like to personally thank you for your Pulsarroll order and it is great you are taking the time to invest in your health and wellbeing.

As a brand we are always trying to improve and give you the best possible customer experience.

We spend countless hours in product development and research putting the needs of you, the customer, at the top of our priority list.

We always appreciate customer feedback and would love to read your reviews and see images using your Pulsarroll across our various social media platforms.

Paul McCabe - Founder and CEO

THE NEXT GENERATION OF VIBRATION TRAINING AND RECOVERY

OFFICIAL SUPPLIERS

WWW.PULSEROLL.COM/WARRANTY

YOU MUST REGISTER YOUR PRODUCT WITHIN 30 DAYS TO ACTIVATE YOUR 1 YEAR WARRANTY

-  Built in rechargeable battery
-  Scientifically proven
-  Reduce next day soreness
-  Increase blood flow
-  Used by Physios & therapists
-  Reduce muscular pain
-  Multiple speed settings
-  Approved by physios and pro athletes
-  Rapid training recovery
-  Flush away lactic acid
-  Portable and lightweight
-  Faster recovery from injury

The Pulseroll Vibrating Products are the Next Generation of warm-up and recovery aids for all levels of athletes and fitness enthusiasts. The combination of pressure and vibration helps to reduce muscle tightness, increase flexibility and reduce next day soreness (DOMS) after your training session.

Bringing pain relief to longstanding muscular pains to areas such as the lower back, legs, arms, feet, neck and shoulders. The vibration creates heat that forces fresh blood into the affected areas which then decreases recovery time from workouts.

Pulseroll works similar to a power plate by gently massaging your sore areas (trigger points) and providing you with the benefits and relaxing feeling that you would receive from a deep tissue sports massage. This will help compliment physio & massage sessions while being able to stay in the comfort of your own home.

**PLEASE CHARGE THE MASSAGE GUN FOR 2 HOURS
BEFORE USING FOR THE FIRST TIME**

OPERATING YOUR PULSEROLL PERCUSSION MASSAGE

- Press the power button once to display how much power the gun has remaining
- Press the power button to turn the gun on and hold for 3 seconds to power off.
- To increase the vibration level of the (Level 1-4), press the power button.
- Each time you press the power button, the vibration level will increase by 1 level

Vibration levels

Level 1 1300 RPM - Low power mode

Level 2 1800 RPM - Medium power mode

Level 3 2500 RPM - High power mode

Level 4 3300 RPM - Ultra high power mode

Included in the box:

Pulserroll Percussion Massage Gun,
user manual,
warranty card,
product charger,
international plug adapters,
Pulserroll Carry Case.

ACCESSORIES INCLUDED

Round Head

For full body muscle massage

Impact Level

Flat Head

For full body muscle massage

Impact Level

Spinal Head

For spinal massage, stimulating massage

Impact Level

Pneumatic Head

For soft tissue and full body muscle massage

Impact Level

Round Head

For full body muscle massage

Impact Level

Bullet Head

For pain point muscle massage

Impact Level

ATTACHING THE MESSAGE HEAD


Fig.1 (Installation)



Fig.2 (Completed Result)

Turn off the machine, apply your chosen attachment in the direction shown in Fig.1. Insert into the groove side of the machine (blue part in Fig.1). Push until connection feels secure. Completed result shown in Fig.2.

Chest:

To relax the chest:

Hold product with reverse grip
move the device slowly along
the muscle in the directions shown

To activate the pectoral muscle:

Use the same action but greatly increase
speed of movement whilst moving up through
power levels.

Duration:

Up to 2 minutes

Attachment:


Small Round

Duration:

Maximum 30 seconds

Attachment:

Small Round

Deltoids

To relax the area:

Hold product, again, with a reverse grip
move the device slowly along
the muscle in the directions shown

To activate the deltoid muscle:

Use the same action but greatly increase
speed of movement whilst moving up through
power levels.

Duration:

Up to 2 minutes

Attachment:


Small Round

Duration:

Maximum 30 seconds

Attachment:

Small Round

Forearms

To relax the area:

With reverse grip move the device slowly along
the muscle in the directions shown.
Lighter pressure advised.

To activate the Forearms:

Use the same action but increase speed of
movement and pressure whilst progressing
through power levels.

Duration:

Up to 2 minutes

Attachment:

Small or Large Round

Duration:

Maximum 30 seconds

Attachment:

Small or Large Round,
Flat head

Hands

To relax the area:

This time using a normal forward facing grip.
Use a lighter pressure whilst moving along the
entire area.

To activate the Hands:

Use the same action throughout the whole area
but increase speed of movement and pressure
whilst progressing through power levels.

Duration:

Up to 2 minutes

Attachment:


Pneumatic Head

Duration:

Maximum 30 seconds

Attachment:

Pneumatic Head


Lower Back

To relax the area:

You may find holding the product towards the base easiest when reaching around the hips.
Move the device slowly along the area in the directions shown

To activate the muscle group:

Use the same method but increasing speed of movement and progressing through power levels.


Upper arm and the Bicep

To relax the area:

We recommend a reverse grip
move the device slowly along the area
in the directions shown

To activate the muscle group:

Use the same action but greatly increase speed of movement whilst moving up through power levels.

Upper Back

To relax the area:

With the base grip move the device slowly along the muscle in the directions shown.
Lighter pressure advised.

To activate the Forearms:

Use the same action but increase speed of movement and pressure whilst progressing through power levels as necessary.

Hamstrings

To relax the area:

Using the reverse facing grip.
Use a lighter pressure whilst moving along the entire area.

To activate the Hands:

Repeat the exercise throughout the whole area but increase speed of movement and pressure whilst progressing through power levels.

Duration:

Up to 2 minutes

Attachment:

Flat Head, spinal head

Duration:

Maximum 30 seconds

Attachment:

Flat Head, spinal head

Duration:

Up to 2 minutes

Attachment:

Small Round

Duration:

Maximum 30 seconds

Attachment:

Small Round

Duration:

Up to 2 minutes

Attachment:

Pneumatic Head

Duration:

Maximum 30 seconds

Attachment:

Pneumatic Head, Flat head

Duration:

Up to 2 minutes

Attachment:


Large round

Duration:

Maximum 30 seconds

Attachment:

Large Round, Flat head

Hip Flexors

To relax the area:

Using a reverse grip facing the product towards the area.

Move the device slowly along the area in the directions shown

To activate the muscle group:

Utilise the same exercise this time increasing speed of movement whilst increasing level to suit.

IT Band

To relax the area:

We recommend a reverse grip with the product upside down move the device slowly along the area in the directions shown

To activate the band:

Use the same action but increase speed of movement whilst cycling through power levels, increasing pressure as necessary.

Quads

To relax the area:

Using the standard reverse grip face the product towards the area. Move slowly along the muscle in the directions shown. Lighter pressure advised.

To activate the muscle group:

Using the same action but increasing speed of movement and pressure whilst progressing through power levels over a shorter duration.

Adductors

To relax the area:

Using the reverse facing grip.

Use a lighter pressure whilst moving along the entire area.

To activate the adductors:

Repeat the exercise throughout the whole area but increase speed of movement and pressure whilst progressing through power levels.

Duration:

Up to 2 minutes

Attachment:

Small round

Duration:

Maximum 30 seconds

Attachment:

Small round

Duration:

Up to 2 minutes

Attachment:

Flat Head

Duration:

Maximum 30 seconds

Attachment:

Flat Head

Duration:

Up to 2 minutes

Attachment:

Large Ball

Duration:

Maximum 30 seconds

Attachment:

Large Ball

Duration:

Up to 2 minutes

Attachment:

Large round

Duration:

Maximum 30 seconds

Attachment:

Large Round

Calves

To relax the area:

Hold the unit at the base to allow for easiest motion in the desired direction.

Cover the area using lighter pressure in the directions shown.

To activate the muscle group:

Utilise the same exercise over a shorter period increasing pressure and speed of movement.

Duration:

Up to 2 minutes

Attachment:

Pneumatic Head

Duration:

Maximum 30 seconds

Attachment:

Pneumatic Head

FEET

To relax the area:

We recommend a reverse grip with the product upside down move the device slowly along the area in the directions shown apply lighter pressure.

To activate the band:

Use the same action but increase speed of movement change Attachment where necessary.

Duration:

Up to 2 minutes

Attachment:

Pneumatic or large Round

Duration:

Maximum 30 seconds

Attachment:

Small round or Bullet

Glutes

To relax the area:

Using the reverse facing grip.

Use a lighter pressure whilst moving along the entire area.

To activate the adductors:

Repeat the exercise throughout the whole area but increase speed of movement and pressure whilst progressing through power levels.

Duration:

Up to 2 minutes

Attachment:

Large round

Duration:

Maximum 30 seconds

Attachment:

Large Round

Rhomboids

To relax the area:

Using the reverse facing grip.

Use a lighter pressure whilst moving along the entire area.

To activate the muscle:

Repeat the exercise throughout area but increase speed of movement and pressure and progress through power levels to suit

Duration:

Up to 2 minutes

Attachment:

Large round


Duration:

Maximum 30 seconds

Attachment:

Large Round

Tennis Elbow

To ease pain in the area:

Hold the product in a reverse grip facing towards the area. Slowly move the product around the entire elbow area along with the bottom of the tricep.


Use a lighter pressure for the Duration shown.

Duration:

10 Seconds Maximum

Attachment:

Pneumatic Head

Shinsplints

To ease pain in the area:

We recommend a reverse grip move the device slowly up and down the area in the directions show using a moderate pressure

Duration:

30 Seconds Maximum

Attachment:

Pneumatic Head

Plantar Fascistic

To ease pain in the area:

Recommending a reverse grip and a circular motion. Treat the foot as 3 sections

1. Ball
2. Arch
3. Heel


Use moderate pressure for the duration Recommended.

Duration:

Up to 20 seconds in each area of the foot

Attachment:

Pneumatic Head

Carpel tunnel

To ease the pain in this area:

Using the reverse facing grip. Use a lighter pressure whilst moving along the entire area forearm to wrist.

Duration:

30 seconds Maximum

Attachment:

Pneumatic Head

**TO SEE HOW TO USE VIDEOS PLEASE
VISIT OUR WEBSITE
WWW.PULSEROLL.COM/HOW-TO-USE**

IMPORTANT SAFETY INSTRUCTIONS

Be sure to read the instructions in this section carefully before use. We also recommend you keep this instruction manual for future reference.

WARNINGS & CAUTION

- Do not place in or drop into water or other liquid.
- Close supervision is necessary when this unit is used by/on children under 10 years old and ailing or disabled persons.
- The product should never be used by people with medical conditions or by people under medical care without first consulting their doctor.
- Check your doctor as well if pregnant, breast implants or pacemakers are present, if you have circulatory disorders, wounds, bruises, inflammations or varicose veins.
- Use the product for its intended use only as described in this manual.
- Always unplug the product after charging or when preparing to use.
- Do not continuously use for more than 60 minutes. Recommended to stop for 20 minutes then restart to ensure LONGER LIFE of the unit.
- While in operation please keep the product away from the neck area.

MAINTENANCE & CLEANING

Unplug before cleaning and then clean only with a towel or soft cloth. Do not use brushes, aggressive cleaning agents, petrol or thinners. Gently dry the unit with a clean dry cloth.

Never put the product in water to clean it and always make sure that the water or any other fluid cannot get into the product.

Do not use the product if it has a damaged cable or plug; if it is not working properly then contact Pulseroll to discuss repair or replacement. The integrated batteries are not designed for exchange or removal by the customer. If problems with the power supply should occur, contact Pulseroll. If signs of water inside the product are present then the warranty may then be void.

BATTERY CHARGING

- Your Pulseroll operates on Lithium Ion batteries (7.4V,2200mAh)
- The product will not charge unless it has been turned off.
- Charge the product before using it for the first time until the battery indicator shows it's fully charged.
- To charge the batteries fully, the charging time will take approx 3 hours. When fully charged you can use for up to 4 hours (level 1) or 2 hours (level 4)
- If the vibration intensity ceases significantly then turn off and fully recharge.
- The product may be recharged at any time regardless of the battery level.
- The product needs to be fully recharged at least once every 3 months.
- Do not ever plug in or remove the power supply plug with wet hands; it may cause electric shock.

POWER SUPPLY

- Use the product only with the voltage specified on the unit. Using a voltage higher than that which is recommended may cause fire or electric shock.
- Only insert the power supply unit into a socket when the product has been switched off.
- Only use the power supply unit supplied for your Pulseroll.
- Charge the product before using it for the first time until the battery indicator shows it's fully charged.
- Do not cut, kink, otherwise damage or modify the power supply cord.
- Keep the power supply unit, cable and product far away from heat, hot surfaces, moisture and liquids.
- Do not touch the power supply unit when you are standing in water and always do take hold of the plug with dry hands.
- Do not touch the product or power supply if it falls into water when plugged into charge. Pull out the power plug and the power supply unit from the socket immediately.

PRODUCT WARRANTY / REPAIR & REPLACEMENT

Pulseroll warrants this product to be free of defects in material or workmanship for a period of one year from the date of original retail purchase. This warranty is available only to the original purchaser, and only if purchased from an authorized Pulseroll distributor or reseller. If your product is found to be defective during the warranty period, Pulseroll will at its discretion, either repair the product or replace the product with a new or factory reconditioned or substantially equivalent model. This warranty applies only to the external shell of the product, the internal components and charger.

THE WARRANTY IS SUBJECT TO BELOW TERMS:

- Pulseroll must be notified of the defect.
- Proof of purchase must be presented.
- Warranty will be void if product is modified, misused or repaired by an unauthorised person.
- Warranty after repair will not be extended beyond the original one year period.
- All replacement parts will be new or reconditioned.
- The customer shall pay the postal charges.
- Warranty must be activated via pulseroll.com/warranty

THE WARRANTY DOES NOT COVER

- Damage to the by improper power supply voltage, accident, fire, floods or acts of nature.
- Damage resulting from misuse and abnormal use or conditions; improper storage. transport and rough handling; exposure to moisture dampness or an environment where corrosion may occur; neglect and abuse; food and any liquid spillage.

PULSEROLL EDUCATION

Pulseroll education is ideal for customers and professionals recommending products to clients, even for someone who has purchased the products and simply would like to know more. This ensures that everyone that takes the education knows the very best ways to get the maximum benefits out of the roller and peanut ball, that will be ultimately passed on to customers.

Alternatively visit pulseroll.com for further tutorials and how to guides.

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