SCIENTIFICALLY
PROVEN

pulseroll

VIBRATING FOAM ROLLER PRO MANUAL

OFFICIAL SUPPLIERS

BRITISH
WEIGHT
LIFTING

BRITISH
FENCING

DEAR CUSTOMER

I would like to personally thank you for your Pulseroll order and it is great you are taking the time to invest in your health and wellbeing.

As a brand we are always trying to improve and give you the best possible customer experience.

We spend countless hours in product development and research putting the needs of you, the customer, at the top of our priority list.

We always appreciate customer feedback and would love to read your reviews and see images using your Pulseroll across our various social media platforms.

Paul McCabe - Founder and CEO

THE NEXT GENERATION OF VIBRATION TRAINING AND RECOVERY

OFFICIAL SUPPLIERS

WWW.PULSEROLL.COM/WARRANTY

YOU MUST REGISTER YOUR PRODUCT WITHIN 30 DAYS TO ACTIVATE YOUR 1 YEAR WARRANTY

-  Built in rechargeable battery
-  Scientifically proven
-  Reduce next day soreness
-  Increase blood flow
-  Used by Physios & therapists
-  Reduce muscular pain
-  Multiple speed settings
-  Approved by physios and pro athletes
-  Rapid training recovery
-  Flush away lactic acid
-  Portable and lightweight
-  Faster recovery from injury

The Pulseroll Vibrating Products are the Next Generation of warm-up and recovery aids for all levels of athletes and fitness enthusiasts. The combination of pressure and vibration helps to reduce muscle tightness, increase flexibility and reduce next day soreness (DOMS) after your training session.

Bringing pain relief to longstanding muscular pains to areas such as the lower back, legs, arms, feet, neck and shoulders. The vibration creates heat that forces fresh blood into the affected areas which then decreases recovery time from workouts.

Pulseroll works similar to a power plate by gently massaging your sore areas (trigger points) and providing you with the benefits and relaxing feeling that you would receive from a deep tissue sports massage. This will help compliment physio & massage sessions while being able to stay in the comfort of your own home.

SAMPLE EXERCISES

Use Pulseroll on your sore muscle areas, hold for 3-5 minutes on the most painful trigger point and then continue to roll. Repeat several times each session and on a daily basis for great long term results and benefits. The intensive Level 4 vibration will help achieve the best results.

**TO SEE HOW TO USE VIDEOS PLEASE
VISIT OUR WEBSITE**

WWW.PULSEROLL.COM/HOW-TO-USE

HAMSTRINGS

In a seated position, place the Pulseroll under the back of the upper leg. You can either cross your other leg over or have both extended over the foam roller. Place your hands by your sides to support your weight. Slowly and deliberately move from the knee to underneath the glute finding trigger points as you go.

LOWER BACK

Place the Pulseroll under your lower back in a seated position raising your hips off the floor. Hold your elbows on the floor and keep your weight onto your lower back.

ILIOTIBIAL TRACT (IT BAND)

Lie on your side and place the Pulseroll under the side of the leg, cross your other leg in front with your knee bent and foot flat to the floor. Place as much weight onto the roller that is tolerable. Slowly move towards the knee working your way through to help massage your trigger points.

CALVES

Place the Pulseroll underneath your lower leg, then cross your other leg over. Place your hands by your sides and press down to raise your hips off the floor, pushing your weight against the calf muscle. Move the leg from side-to-side while keeping the roller static to find tight areas before moving up towards the knee.

SHOULDERS

Place the Pulseroll underneath your shoulder blades and lift your hips off the ground to create more pressure if needed with your hands in a relaxed position by your sides. Keep your knees bent and feet flat to the floor. Move slowly to find tight areas and stick on them until the tightness eases.

GLUTES

Sit on the Pulseroll and shift your weight to one side. slowly move from the top of your hamstring up to your bottom area, finding tight areas and then sticking on them while keeping your hands on the floor for support. Repeat for the other side.

QUADRICEPS

Lying face down to the floor, position the Pulseroll underneath both quads, support your weight using your forearms. Slowly roll from above the knee to below the hips to massage the area.

SINGLE QUADRICEP

Lying face down to the floor, place the Pulseroll underneath the quad. Position the other leg outward to your side. Use your forearms to support your weight. Roll over the quadricep, keeping your forearms firmly to the floor. Repeat for the other leg.

PLANK

Place the Pulseroll underneath the forearms with the feet a shoulder width apart then hold the plank position keeping the core engaged throughout. Keep the head and spine neutral throughout the exercise. Aim for 60 seconds.

PLANK

Place the Pulseroll underneath both feet keeping the core engaged throughout. Stabilise your weight on the elbows making sure they are positioned underneath your shoulder. Aim for 60 seconds.

FEET


Sitting on a chair with your lower back supported put one foot at a time on the Pulseroll. Put your foot down firmly while changing angles allowing the central, inner and outer parts of your foot to massage against the roller. Press your other foot on top for added pressure if required. Repeat with other foot.

**PLEASE CHARGE THE MASSAGE GUN FOR 2 HOURS
BEFORE USING FOR THE FIRST TIME**

OPERATING YOUR PULSEROLL

- Press the power button for a short time to turn on and for 3 seconds to then turn off.
- To increase the vibration level of the (Level 1 - 5), press the power button.
- Each time you press the power button, the vibration level will increase by 1 level.
- Level 4 will produce a vibrating rhythmic pulse level.
- The timer will automatically stop after 10 minutes of constant use.

- 1** Press power once for battery level indicator
- 2** LED lights to indicate the intensity level
- 3** Power On/Off Button
- 4** Fan
- 5** Connection for the Power Supply to recharge

- Level **1** 1500 RPM - Low Power Mode
- Level **2** 2200 RPM - Medium Power Mode
- Level **3** 2700 RPM - Moderate Mode
- Level **4** 1500 - 3200 RPM - Rhythm Mode
- Level **5** 4300 RPM - High Power Mode

FINGER REMOTE FOR VIBRATING PRODUCTS

- 1 Turn ON / OFF
Change vibration speed
- 2 Finger strap
- 3 Battery slot

The Pulseroll finger remote allows you to turn on / off and cycle through the vibrating speeds of your Pulseroll products without the need to adjust position or stretch awkwardly during your recovery or warm up routine

IMPORTANT SAFETY INSTRUCTIONS

Be sure to read the instructions in this section carefully before use. We also recommend you keep this instruction manual for future reference.

WARNINGS & CAUTION

- Do not place in or drop into water or other liquid.
- Close supervision is necessary when this unit is used by/on children under 10 years old and ailing or disabled persons.
- The product should never be used by persons with medical conditions or by persons under medical care without first consulting their doctor.
- Check your doctor as well if pregnant, breast implants or pacemakers are present, if you have circulatory disorders, wounds, bruises, inflammations or varicose veins.
- Use the product for its intended use only as described in this manual.
- Always unplug the product after charging or when preparing to use.
- Do not continuously use for more than 60 minutes. Recommended to stop for 20 minutes then restart to ensure LONGER LIFE of the unit.
- While in operation please keep the product away from the neck area.

MAINTENANCE & CLEANING

Unplug before cleaning and then clean only with a towel or soft cloth. Do not use brushes, aggressive cleaning agents, petrol or thinners. Gently dry the unit with clean dry cloth

Never put the product in water to clean it and always make sure that the water or any other fluid cannot does not get into the product.

Do not use the product if it has a damaged cable or plug; if it is not working properly then contact Pulseroll to discuss repair or replacement. The integrated batteries are not designed for exchange or removal by the customer. If problems with the power supply should occur, contact Pulseroll. If signs of water inside the product are present then the warranty may then be void.

BATTERY CHARGING

- Your Pulseroll operates on Lithium Ion batteries (7.4V, 2200mAh)
- The product will not charge unless it has been turned off.
- Charge the product before using it for the first time until the battery indicator shows it's fully charged.
- To charge the batteries fully, the charging time will take approx. 3 Hours. When fully charged you can use for up to 6 Hours (level 1) or 2 Hours (level 4).
- If the vibration intensity ceases significantly then turn off and fully recharge.
- The product may be recharged at any time regardless of the battery level.
- The product needs to be fully recharged at least once every 3 months.

POWER SUPPLY

- Use the product only with the voltage specified on the unit. Using a voltage higher than that which is recommended may cause fire or electric shock.
- Only insert the power supply unit into a socket when the product has been switched off.
- Only use the power supply unit supplied for your Pulseroll.
- Charge the product before using it for the first time until the battery indicator shows it's fully charged.
- Do not cut, kink, otherwise damage or modify the power supply cord.
- Keep the power supply unit, cable and product far away from heat, hot surfaces, moisture and liquids.
- Do not touch the power supply unit when you are standing in water and always do take hold of the plug with dry hands.
- Do not touch the product or power supply if it falls into water when plugged into charge. Pull out the power plug and the power supply unit from the socket immediately.

PRODUCT WARRANTY / REPAIR & REPLACEMENT

Pulseroll warrants this product to be free of defects in material or workmanship for a period of one year from the date of original retail purchase. This warranty is available only to the original purchaser, and only if purchased from an authorized Pulseroll distributor or reseller. If your product is found to be defective during the warranty period, Pulseroll will at its discretion, either repair the product or replace the product with a new or factory reconditioned or substantially equivalent model. This warranty applies only to the external shell of the product, the internal components and charger.

THE WARRANTY IS SUBJECT TO BELOW TERMS:

- Pulseroll must be notified of the defect.
- Proof of purchase must be presented.
- Warranty will be void if product is modified, misused or repaired by an unauthorised person.
- Warranty after repair will not be extended beyond the original one year period.
- All replacement parts will be new or reconditioned.
- The customer shall pay the postal charges.
- Warranty must be activated via pulseroll.com/warranty

THE WARRANTY DOES NOT COVER

- Damage to the by improper power supply voltage, accident, fire, floods or acts of nature.
- Damage resulting from misuse and abnormal use or conditions; improper storage. transport and rough handling; exposure to moisture dampness or an environment where corrosion may occur; neglect and abuse; food and any liquid spillage.

PULSEROLL EDUCATION

Pulseroll education is ideal for customers and professionals recommending products to clients, even for someone who has purchased the products and simply would like to know more. This ensures that everyone that takes the education knows the very best ways to get the maximum benefits out of the roller and peanut ball, that will be ultimately passed on to customers.

Alternatively visit pulseroll.com for further tutorials and how to guides.

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